

**Canon**

CANON MEDICAL SYSTEMS DO BRASIL



UNIVERSIDADE  
FEDERAL  
DE PERNAMBUCO

DCM Uploader: um sistema de informação para garantia de confidencialidade, integridade e disponibilidade no tráfego de imagens médicas baseadas no padrão DICOM

**Resultados da Sprint 05b**

25/04/2023

- Investigar implementação da comunicação site to site
  - Pydicom
  - Orthanc
  - Dcm4chee
- Estruturar arquitetura
  - 2 MVP

1. Pydicom
2. Dcm4chee
3. Orthanc

# Resultados - Pydicom/Pynetdicom 1/6

```
pydicom-imagem.py - /home/brunodantas/Teste-PyNetDicom/pydicom-imagem.py
File Edit Format Run Options Window Help

import pydicom
from pydicom.data import get_testdata_file

filename = get_testdata_file("OBXXXX1A_rle.dcm")
ds = pydicom.dcmread(filename)
ds.PatientName
print(ds)
```

```
brunodantas@bruno-dantas:~/pydicom-data/data_store$ cd data
brunodantas@bruno-dantas:~/pydicom-data/data_store/data$ dir
693_J2KR.dcm      RG3_J2KR.dcm      US1_J2KI.dcm
693_UNCI.dcm      RG3_UNCI.dcm      US1_J2KR.dcm
693_UNCR.dcm      RG3_UNCR.dcm      US1_UNCI.dcm
JPEG-LL.dcm       SC_rgb.dcm         US1_UNCR.dcm
JPEG2000_UNC.dcm  SC_rgb_16bit.dcm  bad_sequence.dcm
JPGLosslessP14SV1_1s_1f_8b.dcm  SC_rgb_16bit_2frame.dcm  color-pl.dcm
MR-SIEMENS-DICOM-WithOverlays.dcm  SC_rgb_2frame.dcm  color-px.dcm
MR2_J2KI.dcm      SC_rgb_32bit.dcm  color3d_jpeg_baseline.dcm
MR2_J2KR.dcm      SC_rgb_32bit_2frame.dcm  eCT_Supplemental.dcm
MR2_UNCI.dcm      SC_rgb_dcmrk_ebcr_dcmd.dcm  emri_small.dcm
MR2_UNCR.dcm      SC_rgb_dcmrk_ebcyn1_dcmd.dcm  emri_small_RLE.dcm
OBXXXX1A.dcm      SC_rgb_dcmrk_ebcyn2_dcmd.dcm  emri_small_big_endian.dcm
OBXXXX1A_2frame.dcm  SC_rgb_dcmrk_ebcynp_dcmd.dcm  emri_small_jpeg_2k_lossless.dcm
OBXXXX1A_expb.dcm  SC_rgb_dcmrk_ebcys2_dcmd.dcm  emri_small_jpeg_2k_lossless_too_short.dcm
OBXXXX1A_expb_2frame.dcm  SC_rgb_dcmrk_ebcys4_dcmd.dcm  emri_small_jpeg_ls_lossless.dcm
OBXXXX1A_rle.dcm  SC_rgb_expb.dcm  explicit_VR-UN.dcm
OBXXXX1A_rle_2frame.dcm  SC_rgb_expb_16bit.dcm  gdcms-US-ALOKA-16.dcm
OT-PAL-8-face.dcm  SC_rgb_expb_16bit_2frame.dcm  gdcms-US-ALOKA-16_big.dcm
RG1_J2KI.dcm      SC_rgb_expb_2frame.dcm  liver.dcm
RG1_J2KR.dcm      SC_rgb_expb_32bit.dcm  liver_expb.dcm
RG1_UNCI.dcm      SC_rgb_expb_32bit_2frame.dcm  mlut_18.dcm
RG1_UNCR.dcm      SC_rgb_gdcm2k_uncompressed.dcm  vlut_04.dcm
RG3_J2KI.dcm      SC_ybr_full_uncompressed.dcm
brunodantas@bruno-dantas:~/pydicom-data/data_store/data$ ~
```

```
IDLE Shell 3.8.10
File Edit Shell Debug Options Window Help

TypeError: dcmread: Expected a file path or a file-like, but got NoneType
>>>
===== RESTART: /home/brunodantas/Teste-PyNetDicom/pydicom-imagem.py =====
Squeezed text (100 lines).
>>>
===== RESTART: /home/brunodantas/Teste-PyNetDicom/pydicom-imagem.py =====
Dataset.file_meta -----
(0002, 0000) File Meta Information Group Length  UL: 230
(0002, 0001) File Meta Information Version       OB: b'\x00\x01'
(0002, 0002) Media Storage SOP Class UID        UI: Ultrasound Image Storage
(0002, 0003) Media Storage SOP Instance UID     UI: 1.3.46.670589.14.1000.210.2.199999.20110
525185628.1.0
(0002, 0010) Transfer Syntax UID                UI: RLE Lossless
(0002, 0012) Implementation Class UID          UI: 1.2.826.0.1.3680043.2.1143.107.104.103.1
15.2.8.4
(0002, 0013) Implementation Version Name       SH: 'GDCM 2.8.4'
(0002, 0016) Source Application Entity Title    AE: 'gdcmsconv'
-----
(0008, 0005) Specific Character Set             CS: 'ISO_IR 100'
(0008, 0008) Image Type                        CS: ['ORIGINAL', 'PRIMARY', 'OBSTETRICAL']
(0008, 0016) SOP Class UID                     UI: Ultrasound Image Storage
(0008, 0018) SOP Instance UID                  UI: 1.3.46.670589.14.1000.210.2.199999.20110
525185628.1.0
(0008, 0020) Study Date                        DA: '20110525'
(0008, 0022) Acquisition Date                  DA: '20110525'
(0008, 0023) Content Date                      DA: '20110525'
(0008, 002a) Acquisition DateTime              DT: '20110525145628.350000'
(0008, 0030) Study Time                        TM: '142825.000000'
(0008, 0032) Acquisition Time                  TM: '145628.350000'
(0008, 0033) Content Time                      TM: '145628.350000'
(0008, 0050) Accession Number                  SH: ''
(0008, 0060) Modality                          CS: 'US'
(0008, 0070) Manufacturer                     LO: 'Philips Medical Systems'
(0008, 0080) Institution Name                  LO: 'Philips Healthcare'
(0008, 0090) Referring Physician's Name        PN: ''
(0008, 1010) Station Name                      SH: 'OEM-4K7C02TYJWP'
(0008, 1090) Manufacturer's Model Name         LO: 'CX50'
(0010, 0010) Patient's Name                    PN: 'OB^^^^'
(0010, 0020) Patient ID                       LO: '11-05-25-142825'
```

Ln: 239 Col: 4



# Resultados - Pydicom/Pynetdicom 2/6

imagem-dicom.py - /home/brunodantas/pydicom-data/data

File Edit Format Run Options Window Help

```
import numpy as np
import imageio
from pydicom import dcmread
import matplotlib.pyplot as plt
```

```
dicom_ds = dcmread("./0015.DCM")
print(dicom_ds)
```

```
print(dicom_ds.file_meta)
dicom_ds[0x0008,0x0016]
```

```
img_xray = dicom_ds.pixel_array
img_xray.shape
img_xray = dicom_ds.pixel_array
img_xray.shape
```

imagem-dicom.py - /home

File Edit Format Run Options

```
import numpy as np
import imageio
from pydicom import dcmread
import matplotlib.pyplot as plt
```

```
dicom_ds = dcmread("./0015.DCM")
print(dicom_ds)
```

```
print(dicom_ds.file_meta)
dicom_ds[0x0008,0x0016]
```


```
img_xray = dicom_ds.pixel_array
img_xray.shape
img_xray = dicom_ds.pixel_array
img_xray.shape
```

IDLE Shell 3.8.10

File Edit Shell Debug Options Window Help

```
Python 3.8.10 (default, Mar 13 2023, 10:26:41)
[GCC 9.4.0] on linux
Type "help", "copyright", "credits" or "license()" for more information.
>>>
==== RESTART: /home/brunodantas/pydicom-data/data_store/data/imagem-dicom.py ====
Dataset.file_meta -----
(0002, 0000) File Meta Information Group Length    UL: 202
(0002, 0001) File Meta Information Version         OB: b'\x00\x01'
(0002, 0002) Media Storage SOP Class UID          UI: X-Ray Radiofluoroscopic Image Storage
(0002, 0003) Media Storage SOP Instance UID       UI: 1.2.840.113619.2.15.1008000062035011254.825190719.0.31.2.1
(0002, 0010) Transfer Syntax UID                  UI: Explicit VR Little Endian
(0002, 0012) Implementation Class UID             UI: 1.2.840.113619.6.36
(0002, 0013) Implementation Version Name          SH: '1_2_5'
(0002, 0016) Source Application Entity Title       AE: 'ard-demo'
-----
(0008, 0005) Specific Character Set                CS: 'ISO_IR 100'
(0008, 0008) Image Type                            CS: ['ORIGINAL', 'PRIMARY', 'SI
NGLE PLANE']
(0008, 0016) SOP Class UID                          UI: X-Ray Radiofluoroscopic Image Storage
(0008, 0018) SOP Instance UID                      UI: 1.2.840.113619.2.15.1008000062035011254.825190719.0.31.2.1
(0008, 0020) Study Date                            DA: '19960308'
(0008, 0021) Series Date                           DA: '19960308'
(0008, 0022) Acquisition Date                      DA: '19960308'
(0008, 0023) Content Date                          DA: '19960308'
(0008, 0030) Study Time                            TM: ''
(0008, 0032) Acquisition Time                      TM: '105650'
(0008, 0033) Content Time                          TM: '105650'
(0008, 0050) Accession Number                      SH: ''
(0008, 0060) Modality                              CS: 'RF'
(0008, 0070) Manufacturer                          LO: 'GE MEDICAL SYSTEMS'
(0008, 0090) Referring Physician's Name            PN: ''
(0008, 1010) Station Name                          SH: ''
(0008, 1030) Study Description                     LO: '5'
(0008, 103e) Series Description                    LO: ''
(0008, 1050) Performing Physician's Name           PN: '00558747^'
(0008, 1090) Manufacturer's Model Name             LO: 'DRS'
```

Ln: 92 Col: 4



## Resultados - Pydicom/Pynetdicom 3/6

The image shows a terminal window and an IDE (IDLE 3.8.10) running a Python script. The terminal window displays the command `python3 -m pynetdicom echoscp 11112 -v` and its output, including "Accepting Association" and "Association Released". The IDE window shows the script `my_scu.py` and its output, which includes a restart message and the message "Associação estabelecida com a Echo SCP!".

```
brunodantas@bruno-dantas: ~/estudo-ae
brunodantas@bruno-dantas:~/estudo-ae$ idle &
[1] 7957
brunodantas@bruno-dantas:~/estudo-ae$ python -m pynetdicom echoscp 11112 -v
Command 'python' not found, did you mean:
  command 'python3' from deb python3
  command 'python' from deb python-is-python3
brunodantas@bruno-dantas:~/estudo-ae$ python3 -m pynetdicom echoscp 11112 -v
I: Accepting Association
I: Association Released
```

```
my_scu.py - /home/brunodantas/estudo-ae/my_scu.py (3.8.10)
File Edit Format Run Options Window Help
from pynetdicom import AE


ae = AE()
ae.add_requested_context('1.2.840.10008.1.1')
assoc = ae.associate("127.0.0.1", 11112)
if assoc.is_established:
    print('Associação estabelecida com a Echo SCP!')
    assoc.release()
else:
    print('Falha ao associar')
```

```
IDLE Shell 3.8.10
File Edit Shell Debug Options Window Help
Python 3.8.10 (default, Mar 13 2023, 10:26:41)
[GCC 9.4.0] on linux
Type "help", "copyright", "credits" or "license()" for more information.
>>>
===== RESTART: /home/brunodantas/estudo-ae/my_scu.py =====
Associação estabelecida com a Echo SCP!
>>>
```

## Resultados - Pydicom/Pynetdicom 4/6

```
brunodantas@bruno-dantas: X + v
Requesting Association
D: Request Parameters:
D: ===== OUTGOING A-ASSOCIATE-RQ PDU =====
D: Our Implementation Class UID:      1.2.826.0.1.3680043.9.3811.2.0.0
D: Our Implementation Version Name:  PYNETDICOM_200
D: Application Context Name:         1.2.840.10008.3.1.1.1
D: Calling Application Name:         PYNETDICOM
D: Called Application Name:          ANY-SCP
D: Our Max PDU Receive Size:        16382
D: Presentation Context:
D:   Context ID:      1 (Proposed)
D:   Abstract Syntax: =Verification SOP Class
D:   Proposed SCP/SCU Role: Default
D:   Proposed Transfer Syntaxes:
D:     =Implicit VR Little Endian
D:     =Explicit VR Little Endian
D:     =Explicit VR Big Endian
D: Requested Extended Negotiation: None
D: Requested Common Extended Negotiation: None
D: Requested Asynchronous Operations Window Negotiation: None
D: Requested User Identity Negotiation: None
D: ===== END A-ASSOCIATE-RQ PDU =====
D: Accept Parameters:
D: ===== INCOMING A-ASSOCIATE-AC PDU =====
D: Their Implementation Class UID:    1.2.826.0.1.3680043.9.3811.2.0.0
D: Their Implementation Version Name: PYNETDICOM_200
D: Application Context Name:          1.2.840.10008.3.1.1.1
D: Calling Application Name:          PYNETDICOM
D: Called Application Name:           ANY-SCP
D: Their Max PDU Receive Size:        16382
```

•



```
my_scp.py - /home/brunodantas/estudo-ae-segundo/my_scp.py (3.8.10)
File Edit Format Run Options Window Help
from pydicom.uid import ExplicitVRLittleEndian

from pynetdicom import AE, debug_logger
from pynetdicom.sop_class import CTImageStorage

debug_logger()

ae = AE()
ae.add_supported_context(CTImageStorage, ExplicitVRLittleEndian)
ae.start_server(("127.0.0.1", 11112), block=True)
```




## Supported Service Classes

- Application Event Logging Service Class
- Basic Worklist Management Service Class
- Color Palette Query/Retrieve Service Class
- Defined Procedure Protocol Query/Retrieve Service Class
- Display System Management Service Class
- Hanging Protocol Query/Retrieve Service Class
- Implant Template Query/Retrieve Service Class
- Instance Availability Notification Service Class
- Media Creation Management Service Class
- Modality Performed Procedure Step Managemen

1. Pydicom
2. Dcm4chee
3. Orthanc

## Resultados - DCM4chee 1/5 (Sprint 5b)



## POC

REDE => **AE TITLE** = Chave aleatória ou hash de uma senha

- ID único
- Restrição SCP
- Chave privada
- API key
- AE Title Dicom: 16 caracteres ->  $10^{24}$  combinações


## dcm4chee

- STORESCP
- STORESCU (mais rápido que rsync e SSH)

## Resultados - DCM4chee 3/5 (Sprint 5b)

```
root@BROGBA:/home/medscantrack# node medscan hsS24Asd 6000 /canon6000 /db6000 8080
{
  aetitle: 'hsS24Asd',
  scpport: '6000',
  scpfolder: '/canon6000',
  dbfolder: '/db6000',
  apiport: '8080'
}
DB is running embedded
SCP is running at hsS24Asd:6000
API is running on port 8080
```

▼ medscantrack
▼ SCP 6000
GET :8080/node
GET :8080/scpfiles
POST :8080/node
DEL :8080/node





# Resultados - DCM4chee 4/5 (Sprint 5b)

## Teste com dataset 1 - 660 MB

LOCAL - LOCAL: < 3 s

LOCAL - VM-CIN: < 3 min

```
root@BROGBA:/home/medscantrack# node medscan hsS24Asd 6001 /canon6001 /db6001 8081
{
  aetitle: 'hsS24Asd',
  scpport: '6001',
  scpfolder: '/canon6001',
  dbfolder: '/db6001',
  apiport: '8081'
}
DB is running embedded
SCP is running at hsS24Asd:6001
API is running on port 8081
SCU /canon6001/vm-scdedeploy.cin.ufpe.br_6000 --> hsS24Asd@vm-scdedeploy.cin.ufpe.br:6000
done in 02m43.62s
```

LOCAL - AWS: 49 min / 7 min

```
root@BROGBA:/home/medscantrack# node medscan hsS24Asd 6002 /canon6002 /db6002 8082
{
  aetitle: 'hsS24Asd',
  scpport: '6002',
  scpfolder: '/canon6002',
  dbfolder: '/db6002',
  apiport: '8082'
}
DB is running embedded
SCP is running at hsS24Asd:6002
API is running on port 8082
SCU /canon6002/18.231.81.239_80 --> hsS24Asd@18.231.81.239:80
done in 49m25.66s
SCU /canon6002/18.231.81.239_80 --> hsS24Asd@18.231.81.239:80
done in 07m25.22s
```


# Resultados - DCM4chee 5/5 (Sprint 5b)

## Teste com dataset 2 - 5.29 GB

**LOCAL - LOCAL:** < 30 s

**LOCAL - VM-CIN:** 1h30

```
root@BROGBA:/home/medscantrack# node medscan h
{
  aetitle: 'hsS24Asd',
  scpport: '6001',
  scpfolder: '/canon6001',
  dbfolder: '/db6001',
  apiport: '8081'
}
DB is running embedded
SCP is running at hsS24Asd:6001
API is running on port 8081
SCU /canon6001/vm-scdedeploy.cin.ufpe.br_6000
SCU error
done in 20m50.56s
SCU /canon6001/vm-scdedeploy.cin.ufpe.br_6000
SCU error
done in 16m36.23s
SCU /canon6001/vm-scdedeploy.cin.ufpe.br_6000
SCU error
done in 08m25.37s
SCU /canon6001/vm-scdedeploy.cin.ufpe.br_6000
SCU error
done in 11m13.82s
SCU /canon6001/vm-scdedeploy.cin.ufpe.br_6000
SCU error
done in 21m52.09s
SCU /canon6001/vm-scdedeploy.cin.ufpe.br_6000
SCU error
done in 10m23.20s
SCU /canon6001/vm-scdedeploy.cin.ufpe.br_6000
done in 02m0.16s
```



\* Mais rápido que SSH

**LOCAL - AWS:** 2h08


```
root@BROGBA:/home/medscantrack# node medscan hsS24Asd 6002 /ca
{
  aetitle: 'hsS24Asd',
  scpport: '6002',
  scpfolder: '/canon6002',
  dbfolder: '/db6002',
  apiport: '8082'
}
DB is running embedded
SCP is running at hsS24Asd:6002
API is running on port 8082
SCU /canon6002/18.231.81.239_80 --> hsS24Asd@18.231.81.239:80
SCU error
done in 23m31.82s
SCU /canon6002/18.231.81.239_80 --> hsS24Asd@18.231.81.239:80
SCU error
done in 36m40.85s
SCU /canon6002/18.231.81.239_80 --> hsS24Asd@18.231.81.239:80
SCU error
done in 1h08m39.74s
SCU /canon6002/18.231.81.239_80 --> hsS24Asd@18.231.81.239:80
```

Proximos passos

- validação
- gerenciador



1. Pydicom
2. Dcm4chee
3. Orthanc

- Ngrok / Localtunnel



- Orthanc Foo

(<https://drive.google.com/drive/folders/1Dcfx6riwHWYnROytSp5QjLTm7>)







# Resultados - Orthanc 3/4

- Orthanc Bar

(<https://drive.google.com/drive/folders/1Dcfx6riwHWYnROytSp5QjLTm7>)



- Testes concorrentes
  - Estudos pequeno
  - Estudos médio
  - Estudos grande
- Estudo de 5GB \*

- Focar nos diferentes tipos de testes em ambas as soluções
- Comparar os resultados de ambas as soluções
- Melhorar robustez da comunicação em situações adversas
- Simplificar configuração da solução

**Canon**

CANON MEDICAL SYSTEMS DO BRASIL



UNIVERSIDADE  
FEDERAL  
DE PERNAMBUCO

DCM Uploader: um sistema de informação para garantia de confidencialidade, integridade e disponibilidade no tráfego de imagens médicas baseadas no padrão DICOM

**Resultados da Sprint 05a**

18/04/2023

- Analisar diferentes ferramentas de comunicação DICOM
  - Pydicom
  - DCM4che
  - Orthanc
- Elaborar requisitos do sistema
- Propor arquitetura



- Pydicom



Documentação:

- Github: <https://github.com/pydicom/pydicom>
- Biblioteca: <https://pydicom.github.io/pydicom/stable/>
- Versão do Pydicom: 2.3.1

# Resultados - Pydicom

```
brunodantas@bruno-dantas:~/pydicom-data/data_store$ cd data
brunodantas@bruno-dantas:~/pydicom-data/data_store/data$ dir
693_J2KR.dcm          RG3_J2KR.dcm          US1_J2KI.dcm
693_UNCI.dcm          RG3_UNCI.dcm          US1_J2KR.dcm
693_UNCR.dcm          RG3_UNCR.dcm          US1_UNCI.dcm
JPEG-LL.dcm           SC_rgb.dcm            US1_UNCR.dcm
JPEG2000_UNC.dcm      SC_rgb_16bit.dcm      bad_sequence.dcm
JPGLosslessP14SV1_1s_1f_8b.dcm SC_rgb_16bit_2frame.dcm color-pl.dcm
MR-SIEMENS-DICOM-WithOverlays.dcm SC_rgb_2frame.dcm color-px.dcm
MR2_J2KI.dcm          SC_rgb_32bit.dcm      color3d_jpeg_baseline.dcm
MR2_J2KR.dcm          SC_rgb_32bit_2frame.dcm eCT_Supplemental.dcm
MR2_UNCI.dcm          SC_rgb_dcmTk_ebcrc_dcmd.dcm emri_small.dcm
MR2_UNCR.dcm          SC_rgb_dcmTk_ebcyn1_dcmd.dcm emri_small_RLE.dcm
OBXXXX1A.dcm          SC_rgb_dcmTk_ebcyn2_dcmd.dcm emri_small_big_endian.dcm
OBXXXX1A_2frame.dcm   SC_rgb_dcmTk_ebcynp_dcmd.dcm emri_small_jpeg_2k_lossless.dcm
OBXXXX1A_expb.dcm     SC_rgb_dcmTk_ebcys2_dcmd.dcm emri_small_jpeg_2k_lossless_too_short.dcm
OBXXXX1A_expb_2frame.dcm SC_rgb_dcmTk_ebcys4_dcmd.dcm emri_small_jpeg_ls_lossless.dcm
OBXXXX1A_rle.dcm      SC_rgb_expb.dcm       explicit_VR_UN.dcm
OBXXXX1A_rle_2frame.dcm SC_rgb_expb_16bit.dcm gdcm-US-ALOKA-16.dcm
OT-PAL-8-face.dcm     SC_rgb_expb_16bit_2frame.dcm gdcm-US-ALOKA-16_big.dcm
RG1_J2KI.dcm          SC_rgb_expb_2frame.dcm liver.dcm
RG1_J2KR.dcm          SC_rgb_expb_32bit.dcm liver_expb.dcm
RG1_UNCI.dcm          SC_rgb_expb_32bit_2frame.dcm mlut_18.dcm
RG1_UNCR.dcm          SC_rgb_gdcm2k_uncompressed.dcm vlut_04.dcm
RG3_J2KI.dcm          SC_ybr_full_uncompressed.dcm
brunodantas@bruno-dantas:~/pydicom-data/data_store/data$ ~
```

```
Dataset.file_meta -----
(0002, 0000) File Meta Information Group Length  UL: 192
(0002, 0001) File Meta Information Version       OB: b'\x00\x01'
(0002, 0002) Media Storage SOP Class UID         UI: Computed Radiography Image
Storage
(0002, 0003) Media Storage SOP Instance UID       UI: 1.3.6.1.4.1.5962.1.1.9.1.2.
20040826185059.5457
(0002, 0010) Transfer Syntax UID                 UI: JPEG 2000 Image Compression
(Lossless Only)
(0002, 0012) Implementation Class UID            UI: 1.3.6.1.4.1.5962.2
(0002, 0013) Implementation Version Name         SH: 'DCT00L100'
(0002, 0016) Source Application Entity Title      AE: 'CLUNIE1'
-----
(0008, 0005) Specific Character Set               CS: 'ISO_IR 100'
(0008, 0008) Image Type                           CS: ['DERIVED', 'PRIMARY']
(0008, 0012) Instance Creation Date              DA: '20040826'
(0008, 0013) Instance Creation Time              TM: '185542'
(0008, 0014) Instance Creator UID                UI: 1.3.6.1.4.1.5962.3
```

```
pydicom-imagem.py - /home/brunodantas/Teste-Py
File Edit Format Run Options Window Help

import pydicom
from pydicom.data import get_testdata_file
filename = get_testdata_file("693_J2KR.dcm")
filename2 = get_testdata_file("RG1_J2KR.dcm")
ds2 = pydicom.dcmread(filename2)
ds = pydicom.dcmread(filename)
ds.PatientName
ds2.PatientName
print(ds)
print(ds2)
```

```
-----
(0010, 0000) Group Length                        UL: 56
(0010, 0010) Patient's Name                      PN: 'CQ500-CT-310'
(0010, 0020) Patient ID                          LO: 'CQ500-CT-310'
(0010, 0030) Patient's Birth Date                DA: ''
(0010, 0040) Patient's Sex                       CS: ''
(0012, 0062) Patient Identity Removed            CS: 'YES'
(0018, 0000) Group Length                        UL: 360
(0018, 0022) Scan Options                        CS: 'AXIAL MODE'
(0018, 0050) Slice Thickness                     DS: '5.0'
(0018, 0060) KVP                                  DS: '140.0'
(0018, 0088) Spacing Between Slices              DS: '20.0'
(0018, 0090) Data Collection Diameter             DS: '320.0'
(0018, 1020) Software Versions                   LO: 'coreload.81'
(0018, 1100) Reconstruction Diameter             DS: '245.0'
(0018, 1110) Distance Source to Detector         DS: '949.147'
(0018, 1111) Distance Source to Patient          DS: '541.0'
```

# Resultados - Pydicom

imagem-dicom.py - /home/brunodantas/pydicom-data/data

File Edit Format Run Options Window Help

```
import numpy as np
import imageio
from pydicom import dcmread
import matplotlib.pyplot as plt

dicom_ds = dcmread("./0015.DCM")
print(dicom_ds)

print(dicom_ds.file_meta)
dicom_ds[0x0008,0x0016]

img_xray = dicom_ds.pixel_array
img_xray.shape
img_xray = dicom_ds.pixel_array
img_xray.shape
```

imagem-dicom.py - /home

File Edit Format Run Options

```
import numpy as np
import imageio
from pydicom import dcmread
import matplotlib.pyplot as plt
```

```
dicom_ds = dcmread("./0015.DCM")
print(dicom_ds)
```

```
print(dicom_ds.file_meta)
dicom_ds[0x0008,0x0016]
```

```
img_xray = dicom_ds.pixel_array
img_xray.shape
img_xray = dicom_ds.pixel_array
img_xray.shape
```

IDLE Shell 3.8.10

File Edit Shell Debug Options Window Help

Python 3.8.10 (default, Mar 13 2023, 10:26:41)

[GCC 9.4.0] on linux

Type "help", "copyright", "credits" or "license()" for more information.

>>>

==== RESTART: /home/brunodantas/pydicom-data/data\_store/data/imagem-dicom.py ===


Dataset.file\_meta -----

(0002, 0000) File Meta Information Group Length	UL: 202
(0002, 0001) File Meta Information Version	OB: b'\x00\x01'
(0002, 0002) Media Storage SOP Class UID	UI: X-Ray Radiofluoroscopic Image Storage
(0002, 0003) Media Storage SOP Instance UID	UI: 1.2.840.113619.2.15.1008000062035011254.825190719.0.31.2.1
(0002, 0010) Transfer Syntax UID	UI: Explicit VR Little Endian
(0002, 0012) Implementation Class UID	UI: 1.2.840.113619.6.36
(0002, 0013) Implementation Version Name	SH: '1_2_5'
(0002, 0016) Source Application Entity Title	AE: 'ard-demo'

-----

(0008, 0005) Specific Character Set	CS: 'ISO_IR 100'
(0008, 0008) Image Type	CS: ['ORIGINAL', 'PRIMARY', 'SI
(0008, 0016) SOP Class UID	UI: X-Ray Radiofluoroscopic Image Storage
(0008, 0018) SOP Instance UID	UI: 1.2.840.113619.2.15.1008000062035011254.825190719.0.31.2.1
(0008, 0020) Study Date	DA: '19960308'
(0008, 0021) Series Date	DA: '19960308'
(0008, 0022) Acquisition Date	DA: '19960308'
(0008, 0023) Content Date	DA: '19960308'
(0008, 0030) Study Time	TM: ''
(0008, 0032) Acquisition Time	TM: '105650'
(0008, 0033) Content Time	TM: '105650'
(0008, 0050) Accession Number	SH: ''
(0008, 0060) Modality	CS: 'RF'
(0008, 0070) Manufacturer	LO: 'GE MEDICAL SYSTEMS'
(0008, 0090) Referring Physician's Name	PN: ''
(0008, 1010) Station Name	SH: ''
(0008, 1030) Study Description	LO: '5'
(0008, 103e) Series Description	LO: ''
(0008, 1050) Performing Physician's Name	PN: '00558747^'
(0008, 1090) Manufacturer's Model Name	LO: 'DRS'

Ln: 92 Col: 4





# Resultados - Pydicom

The screenshot shows the Visual Studio Code interface. The top menu bar includes File, Edit, Selection, View, Go, Run, Terminal, and Help. The title bar indicates the active file is 'get\_study\_instance\_uid.py' in the 'uploader2' workspace.

The left sidebar contains a SEARCH panel with the following elements:


- SEARCH input field containing 'pydicom'.
- Replace input field.
- A red box highlights the search results, showing 6 results in 5 files. The results are:
  - Dockerfile uploader\_t (1 result): RUN pip3 install pydicom
  - Dockerfile.prod uploader\_t (1 result): RUN pip3 install pydicom
  - get\_ae\_title.py uploader\_t\app\lib\uploader\_t\core.py (2 results): import pydicom, dcm = pydicom.dcmread(filename)
  - get\_study\_description.py uploader\_t\app\lib\uploader\_t\core.py (2 results): import pydicom, dcm = pydicom.dcmread(filename)
  - get\_study\_instance\_uid.py uploader\_t\app\lib\uploader\_t\core.py (2 results): import pydicom, dcm = pydicom.dcmread(filename)

The main editor area displays the code in 'get\_study\_instance\_uid.py'. The code is as follows:

```
1 import pydicom
2 import sys
3 import os
4
5 # Get the file name from the command line
6 # We will extract the AE title from the DICOM file
7 filename = sys.argv[1]
8
9 # Read the file
10 dcm = pydicom.dcmread(filename)
11
12 # Get the metadata from the file and extract the AE_title
13 study_instance_uid = dcm.StudyInstanceUID
14
15 print(study_instance_uid, end="")
16
```

In the code, the 'pydicom' import on line 1 and the 'pydicom.dcmread' call on line 10 are highlighted with a red box.

- Pynetdicom



SCP - Service Class Providers

Mover, Armazenar, Encontra...

SCU - Service Class Users

pynetdicom3

DICOM Upper Layer Protocol for TCP/IP



- dcm4chee
  - Popular framework dicom em java
  - Todas as funcionalidades de um PACS e mais...
  - autenticação, conversão, compressão, criptografia, priorização, redirecionamento, armazenamento, REST, utilities

## Resultados - dcm4chee

- storescp [options] -b [<aet>[@<ip>]:]<port>
- storescu [options] -c <aet>@<host>:<port> [<file>..][<directory>..]


```
root@BROGBA:/home/dcm4chee/bin# sudo ./storescp --accept-unknown -b canon@localhost:6000 --directory /canon
sudo: /etc/sudoers.d is world writable
11:26:10.991 INFO - Start TCP Listener on localhost/127.0.0.1:6000
```

```
.11:32:52.009 INFO - STORESCU->canon(1) << A-RELEASE-RQ
11:32:52.009 DEBUG - STORESCU->canon(1): enter state: Sta7 - Awaiting A-RELEASE-RP PDU
11:32:52.011 INFO - STORESCU->canon(1) >> A-RELEASE-RP
11:32:52.011 INFO - STORESCU->canon(1): close Socket[addr=localhost/127.0.0.1,port=6000,localport=43221]
11:32:52.012 DEBUG - STORESCU->canon(1): enter state: Sta1 - Idle
Sent 643 objects (=316.772MB) in 4.583s (=69.119MB/s)
```

rsync -avhze 'ssh -c aes256-gcm@openssh.com' --compress-level=9 /canon/ /canon2/


```
1.3.6.1.4.1.32839.8440140800298300313739312978940134402
1.3.6.1.4.1.32839.87508610390277843544749621417686490586
1.3.6.1.4.1.32839.90843543623974951824989775779611208119
1.3.6.1.4.1.32839.9665074559205477202144614696690188190
novao
novofile.nada

sent 55.11M bytes  received 5.68K bytes  5.80M bytes/sec
total size is 191.03M  speedup is 3.47
```





- Resultados em 3 ambientes
  - Local
  - VM Cln
  - AWS

## Resultados - Orthanc





# Resultados - Orthanc

- Local



# Resultados - Orthanc



- VM CIn






# Resultados - Orthanc

- AWS



# Resultados - Orthanc



- Investigar implementação da comunicação site to site
- Estruturar arquitetura